



ZHUHAI LITUO BIOTECHNOLOGY CO., LTD.

35TH, YONGAN THREE ROAD, HONGQI TOWN, JINWAN
DISTRICT, ZHUHAI, GUANGDONG, CHINA

STATEMENT

We, Zhuhai Lituo Biotechnology Co., Ltd, hereby specially declare that:

The COVID-19 Antigen Detection Kit (Colloidal Gold) produced by our company detects N protein of SARS-CoV-2, while the mutation of the virus found in England (B.1.1.7 strain), South Africa (B.1.351 strain, B.1.1.529 strain), Brail (P.1 strain) and India (B.1.617 strain, Delta variant B.1.617.2) are mainly on S protein. Therefore, this detection kit still effectively detect the antigen of SARS-CoV-2.

Zhuhai Lituo Biotechnology Co., Ltd. reserves all the right to the interpretation of this letter.



Zhuhai Lituo Biotechnology Co., Ltd.

November 29, 2021